

**A green-extended sensitive , 10-stage, 39mm (1.5") tube**

<b>Applications :</b>	Scintillation counting with green-emitting scintillators or wavelength shifters such as BBQ, CsI(Tl) or BGO.		
<b>Description :</b>	Window :	Material :	lime glass
	Photocathode :		green extended bi-alkali
	Refr. index at 400 nm :		1.54
	Multiplier :	Structure :	linear focused
		Nb of stages :	10
	Mass :		60 g

**Photocathode characteristics**

Spectral range :		290-690	nm
Maximum sensitivity at :		440	nm
Sensitivity ① :			
<input checked="" type="checkbox"/> Luminous :	min.:	90	typ.: 115 $\mu\text{A}/\text{lm}$
Blue :			typ.: 11.5 $\mu\text{A}/\text{lmF}$
Radiant, at 440 nm :			typ.: 90 mA/W

**Characteristics with voltage divider A**

Gain slope (vs supp. volt., log/log) :		7.5	
For an anode blue sensitivity of :		7.5	A/lmF
<input checked="" type="checkbox"/> Supply voltage :	max.:	1250	typ.: 1100 V
	min.:	900	
Gain :		$6.5 \times 10^5$	
<input checked="" type="checkbox"/> Anode dark current ② :	max:	30	typ.: 5 nA
Pulse height resolution <sup>137</sup> Cs ③ :			typ.: 7.2 %
Pulse height resolution <sup>55</sup> Fe ④ :			typ.: 39 %
Peak to valley ratio for <sup>55</sup> Fe :			typ.: 30
Mean anode sensitivity deviation ⑤ :			
long term (16 h) :		1	%
after change of count rate :		1	%
vs temperature between 0 and +40 °C at 420 nm :		+/- 0.1	%/K
Anode sensitivity deviation for a magnetic field of 0.05 mT :		10	%

**Characteristics with voltage divider ⑥ :**

	<b>B</b>	<b>A</b>	
For a supply voltage of :	1500	1500	V
Gain :	$3 \times 10^5$	$6.7 \times 10^6$	
Linearity (2%) of an. current up to :	150	70	mA
Anode pulse ⑦ :			
Rise time :	2.8	2.8	ns
Duration at half height :	7	7	ns
Transit Time :	29	34	ns
Capacitance anode to all :	5		pF

**Recommended voltage divider**

**Type A** for maximum gain

K	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	A	
2	1	1	1	1	1	1	1	1	1	1	1	(total :12)

**Type B** for best timing / linearity compromise

K	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	A	
2	1	1	1	1.25	1.25	1.5	2.25	2.25	2.5	2.25		(total :18.25)



K: photocathode      Dn: dynode      A: anode

**Limiting values**

Anode luminous sensitivity :		max.:	100	A/lmF		
Supply voltage :		max.:	1600	V		
Continuous anode current :		max.:	0.2	mA		
Voltage between :						
	D1 and photocathode :	min.:	100	max.:	500	V
	consecutive dynodes :			max.:	300	V
	anode and D10 :	min.:	40	max.:	300	V
Ambient temperature :						
	short operation (< 30 mn) :	min.:	-30	max.:	+80	°C
	continuous operation & storage :	min.:	-30	max.:	+50	°C

**Notes :**  Characteristic measured and mentioned on the test ticket of each tube.


- ① Luminous sensitivity is measured with a tungsten filament lamp with a colour temperature of  $2856 \pm 5$  K. The blue sensitivity, expressed in A/lmF ("F" as in Filtered) is measured with a tungsten filament lamp with a colour temperature of  $2856 \pm 5$  K. Light is transmitted through an interference filter.
- ② Dark current is measured at ambient temperature, after the tube has been in darkness for approximately 1 min. Lower value can be obtained after a longer stabilisation period in darkness (approx. 30 min.).
- ③ Pulse amplitude resolution for  $^{137}\text{Cs}$  and  $^{57}\text{Co}$  is measured with NaI(Tl) cylindrical scintillator with a diameter of 32 mm and a height of 25 mm. The count rate used is  $\sim 1.0 \times 10^4$  c/s.
- ④ Pulse amplitude resolution for  $^{55}\text{Fe}$  is measured with a NaI(Tl) cylindrical scintillator with a diameter of 25 mm and a height of 1 mm provided with a beryllium window. The count rate used is  $2 \times 10^3$  cps.
- ⑤ The mean pulse amplitude deviation is measured by coupling a NaI(Tl) scintillator to the window of the tube. Long term (16h) deviation is measured by placing a  $^{137}\text{Cs}$  source at a distance from the scintillator such that the count rate is  $\sim 10^4$  c/s, corresponding to an anode current of  $\sim 300$  nA. The mean pulse amplitude deviation after change of count rate is measured with a  $^{137}\text{Cs}$  source at a distance from the scintillator such that the count rate can be changed from  $10^4$  to  $10^3$  c/s, corresponding to an anode current of  $\sim 1$   $\mu\text{A}$  and  $0.1$   $\mu\text{A}$  respectively. Both tests are carried out according to ANSI-N42-9-1972 of IEEE recommendations.
- ⑥ To obtain a peak pulse current greater than that obtainable with divider A, it is necessary to increase the inter-dynode voltage progressively. Divider circuit C is an example of a progressive divider, giving a compromise between gain, speed and linearity. other dividers can be conceived to achieve other compromises. It is generally recommended that the voltage ratio between two successive stages is less than 2.
- ⑦ Measured with a pulse light source, with a pulse duration (FWHM) of approximately 1 ns., the cathode being completely illuminated. The rise time is determined between 10 % and 90 % of the anode pulse amplitude. The signal transit time is measured between the instant at which the illuminating pulse of the cathode becomes maximum, and the instant at which the anode pulse reaches its maximum. Rise time, pulse duration and transit time vary with respect to high tension supply voltage Vht as  $(Vht)^{-1/2}$ .




ref.: 14000030  
sp: short pin  
n: plane of symmetry of the multiplier

K: cathode      Dn: dynode  
A: anode

Typical gain curve



Typical spectral characteristics



Accessories

Socket : FE1112  
Mu-metal shield: MS170