

TOF Detector Performance

- TOF detector specifications are a major factor in tool performance. Key criteria are:
- \blacksquare Timing: faster detector and electronics \rightarrow Improved Mass Resolution
- Output range: high saturation values, better linearity -> Wider dynamic range, single gain working point
- \blacksquare <u>Gain</u>: \rightarrow unambiguous detection of single ion over background noise
- Stability/Lifetime: Longer lifetime \rightarrow Reduced cost of ownership, more stable and simple operation
- El-Mul offers a range of solutions by integrating its ScintiFast[™] technology, fast light sensors and proprietary fast amplifiers into the patented MTOF™ detector concept



MTOF[™] (Magnetic Time of Flight) Concept

- Ions converted to Secondary Electrons (SE) using conversion plate
- Electrons directed to ScintiFastTM scintillator by ExB fields
- Sensor collects photons from scintillator and provides analog signal
- HPD or fast PMT light sensor allows for sub-ns time resolution
- Switchable polarity
- Proprietary fast amplifiers (if needed)
- Nearly 100% detection efficiency
- Flexibility in design to meet space constraints



See also Thermo Fisher Scientific's poster "Characterization of a detection system with high sensitivity and dynamic range for a novel HRAM mass spectrometer" by Johannes Petzoldt et.al..

HPD PMT 6-dynodes PMT 10-dynode the output pulse

- calibration algorithms

MTOF with HPD: Pulse Width and Height

Scintillator Selection

- ScintiFast[™] scintillators can be produced with a range of response times, but there is a tradeoff between decay times and photon yield
- Each TOF system may require a different set of response time and photon yield
- Selection of ScintiFast properties may affect dynamic range
- Selection of ScintiFastTM properties does not affect short- and long-term stability



Versatile, Long-Lifetime, Wide Dynamic Range Detector for TOF Applications

Semyon Shofman, Amit Weingarten, Jonathan Garel, Sasha Kadyshevitch **El-Mul Technologies, Rehovot, Israel**



- Response time / Mass resolution
- Dynamic range
- Gain / noise floor
- Available space
- Hybrid Photo Diode (HPD) and Photo-Multiplier Tubes (PMTs) are compact, low-noise, vacuum-sealed sensors that are suitable for use in MTOF[™] detectors



- Measurements performed on HPD using 80-ps pulse laser
- Linearity is output-limited
- Pulse output linear up to ~200 mV, then monotonic sublinear without pulse distortion up to ~4V
- No dependence on gain
- Overall dynamic range of ~3 orders of magnitude with non-linearity corrections

Sensor Performance

	FWHM (ns)	Max Gain	Operating Voltage
	0.6	1.2E5	8kV
5	0.8	2.5E5	1kV
S	1.1	5E6	1kV

Dynamic range is determined by linearity of 6-dynode PMT are linear up to >6V

10-dvnode PMT range can be extended by

using modified dynode voltage arrangement HPD has the lower range, regardless of gain Actual dynamic range can be extended using



Traces from MTOF detector with HPD including pre-Amplifier operated at maximal gain

FWHM = 0.63 ns (averaged over of 1000 single pulses) Pulse can be shaped by optimizing output impedance Smooth decay without undershoot or significant "ringing" can be achieved without compromising mass resolution

Pulse Height Distribution measured with single electrons on subassembly (ScintiFast + Light-guide + HPD)

- Single electron pulse height = 7.4 ± 2.3 mV
- Ion signal is typically 2-4 electrons = 15-30mV Reliable Single Ion detection above noise level
- Allows to work at lower gain and increase dynamic range



MTOF with PMT: Pulse Width and Height

- Traces from MTOF detector with two PMT types
- FWHM between 1ns for 6-dynodes PMT to 1.8 ns for 10dynodes PMT
- Smooth decay without undershoot or significant "ringing" can be achieved by electronics optimization but pulse is slightly broadened
- Pulse Height Distribution measured with single electrons on subassembly (ScintiFast + Light-guide + PMT)
- Single electron pulse for 6-dynode PMT= 12.4 ± 5.3 mV
- Ion signal is typically 2-4 electrons = 25-50mV
- As with HPD, it allows reliable Single Ion detection above noise level or to decrease gain and increase dynamic range

MTOF with PMT: Linearity and Dynamic Range



El-Mul Technologies Proprietary Information

MTOF with HPD: Linearity and Dynamic Range



MTOF Detectors Configurations

Few examples of TOF detectors based on MTOF technology are shown below

• For each detector, selection of scintillator and sensor, together with the limitations of space and cost were taken into account in the design

Detector	Scintillator	Sensor	Pulse width	Dynamic range
1	Fast	HPD	<0.7ns	Up to 100
2	Semi-fast	6-dynodes PMT	1	2,000
3	Semi-fast	10-dynodes PMT	1.7	2,000
4	"Standard"	10-dynodes PMT	2	>10,000



Multi-ion Timing

- Ion time-spread is mostly due to electric field nonuniformity between entrance grid and convertor, and to convertor non-planarity
- Spread will be larger for larger ion-beam dimensions
- In the upper example, of fast detector (fast scintillator HPD) ions are uniformly distributed over active area of 15*25 mm (worst case scenario)
- Calculation for 10,000 ions with M = 3 kDa and energy 20keV show that FWHM ~ 0.10 ns, and ~90% of positive ions arrive within 0.15 ns, ~100% arrive within 0.25 ns
- The lower figure, an experiment with a slower detector (10-dynode PMT), shows ~10% broadening above 1k
- For several detectors over the parameters range shown, multi-ion signal broadening is typically 10-20%





Light Splitting for Enhanced Dynamic Range

- Modification of the optical system of the MTOF detector enables a wider dynamic range, using a single scintillator and splitting the photons into two channels, where the weak channel collects only a few % of the signal:
- Use the weak channel to collect signals with a large number of ions Use the strong channel to collect signals with a small number of ions

Strong channel: 5ph-e/e \rightarrow

- Each channel includes its own light sensor (HPD or PMT) and electronics
- Ratio between channels can be determined by design dynamic range can be chosen
- Demo unit was scanned with electron beam showing uniform collection in both channels from all scintillator area

Weak channel collects 3%





"strong" channel

64 128 256 512 1024 2048 4096 8192 16384

Output amplitude on 50 Ohm (mV)

-

--- PMT 700V

-0-PMT 900V

- PMT 1000V

Stability and Lifetime

- Full detector, with semi-fast scintillator, and 10 dynodes PMT was measured on El-Mul Benchtop TOF system, PFTBA gaseous source
- Input current 2×10^7 ion/sec $\rightarrow 3$ pA $\rightarrow 0.01 \mu$ C/H
- Detector operated continuously over 1000H showing no significant degradation
- Accumulated to 16 μC input charge (equivalent to ~16 C output charge in "standard" detector)
- Scintillator degradation was also tested using ebeam up to 350 μ C input current and seen to be <10%.





Summary

- Variety of TOF detectors based on El-Mul's patented MTOF [™] concept, incorporating El-Mul's fast ScintiFast™ scintillators are presented
- Selection of the scintillator and sensor types enables detector design for specific performance over wide range of parameters: mass resolution, dynamic range, beam dimensions, available space, cost
- FHWM between <0.7 to 2 ns for single ion pulse was demonstrated</p>
- Broadening for multi-ion signal is typically 10-20% and depends on system parameters
- Dynamic range between 10² of 10⁴ was demonstrated and can be further enhanced by using split optical channels
- Stability and long-lifetime was proven for both fast, narrow dynamic range and for slower, wide dynamic range detectors
- This ultrafast, high gain, wide dynamic range and long lifetime detector represents a significant performance improvement over traditional MCP-based and electron multiplierbased high resolution TOF detectors

Detector Stability



